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

Network Settings

Interface	10Mbps Ethernet	
Current Settings (Static)	IP Address	192.168.33.109
	Subnet Mask	255.255.255.0
	Default Gateway	192.168.33.251
	DNS Server Address	
	DNS Domain Suffix	
Stored Settings	<input checked="" type="radio"/> Obtain an IP address automatically	
	<input type="radio"/> Use the following IP address:	
	IP Address	<input type="text" value="192.168.33.109"/>
	Subnet Mask	<input type="text" value="255.255.255.0"/>
	Default Gateway	<input type="text" value="192.168.33.251"/>
	DNS Server Address	<input type="text"/>
	DNS Domain Suffix	<input type="text"/>

Note: Changes only take effect after [reboot](#).

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Network Settings

Interface	10Mbps Ethernet	
Current Settings (Static)	IP Address	192.168.33.109
	Subnet Mask	255.255.255.0
	Default Gateway	192.168.33.251
	DNS Server Address	
	DNS Domain Suffix	
Stored Settings	<input type="radio"/> Obtain an IP address automatically	
	<input checked="" type="radio"/> Use the following IP address:	
	IP Address	<input type="text" value="192.168.33.109"/>
	Subnet Mask	<input type="text" value="255.255.255.0"/>
	Default Gateway	<input type="text" value="192.168.33.251"/>
	DNS Server Address	<input type="text"/>
	DNS Domain Suffix	<input type="text"/>


Note: Changes only take effect after [reboot](#).

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Stored Settings	<input type="radio"/> Obtain an IP address automatically	
	<input checked="" type="radio"/> Use the following IP address:	
	IP Address	192.168.33.45
	Subnet Mask	255.255.255.0
	Default Gateway	192.168.33.251
	DNS Server Address	
DNS Domain Suffix		

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The image shows a 'Network Settings' dialog box for the '10Mbps Ethernet' interface. It is divided into 'Current Settings (Static)' and 'Stored Settings'. The 'Stored Settings' section shows the 'Use the following IP address' option selected, with the IP address field containing '192.168.33.45'. A red box highlights the 'Submit' button at the bottom left, with a red arrow pointing to it from a separate 'Submit' button callout. The text 'クリック' (Click) is written in red next to the arrow. Below the dialog box, there are 'Submit' and 'Clear Changes' buttons, and a note: 'Note: Changes only take effect after [reboot](#).'

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The image shows an administrator login dialog box titled '192.168.33.109 に接続'. It has a 'Administrator' label and a red text prompt 'root と入力' (Enter 'root'). The 'ユーザー名(U):' (Username) field contains 'root'. The 'パスワード(P):' (Password) field is empty. There is a checkbox for 'パスワードを記憶する(R)' (Remember password) which is unchecked. At the bottom, the 'OK' button is highlighted with a red box, and the text 'クリック' (Click) is written in red next to it. A 'キャンセル' (Cancel) button is also present.

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Contact [Fargo Electronics Inc.](#) for product information, technical support, and downloads.